
# NANOREPRO SARS-CoV-2 Saliva Antigen Test

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NanoRepro AG  
Untergasse 8  
D-35037 Marburg  
www.nanorepro.com

# NANOREPRO SARS-CoV-2 Saliva Antigen Test

Rapid test for the detection of SARS-CoV-2 Antigen

## INSTRUCTION FOR USE

### INTRODUCTION

In December 2019, the novel respiratory disease (COVID-19) caused by the coronavirus (SARS-CoV-2) was reported in Wuhan, China.<sup>1,2</sup> According to WHO, most of the people infected with SARS-CoV-2 have mild to moderate respiratory diseases, fever, cough and recover without treatment. However, people with weak immune systems, such as the elderly or people with previous illnesses (e.g., cardiovascular disease, diabetes, chronic respiratory diseases, cancer, etc.) are more likely to develop a serious illness that can lead to the death of the infected person.<sup>3</sup>

The NanoRepro SARS-CoV-2 Saliva Antigen Test is intended for the qualitative detection of SARS-CoV-2 antigen in human saliva samples from novel coronavirus suspected person with symptoms within 7 days from onset. The test is only for professional use. Positive result of the antigen test can be used for early isolation of patients with suspected infection, but it cannot be used as diagnosis basis of SARS-CoV-2 infection. Negative results do not rule out SARS-CoV-2 infection and basis for treatment. Further nucleic acid detection should be carried out for suspected population whose antigen test result is positive or negative.

The NanoRepro SARS-CoV-2 Saliva Antigen Test is an immunochromatography assay which detects SARS-CoV-2 nucleocapsid antigen in the samples with the help of the double antibody sandwich method. If there is virus antigen presence in the sample, it binds with the corresponding colloidal gold antibody. This complex "migrates" across the membrane and binds to the monoclonal antibody at the Test line (T). This creates a visible red line, which indicates a positive result. However, if the sample does not contain any antigen, then the complex cannot be formed and thus no reddish line forms in the Test line (T). Regardless of whether the sample contains antigen or not, a reddish line forms in the Control line (C).

An explanation of how to read and interpret the test results requires a thorough reading of the instructions for use. It is important to fully understand the entire instruction pamphlet before performing the test.

### TEST CONTENTS

- 10 test cassettes with 1 pipette in a sealed pouch (SARS-CoV-2 Ag)
- 10 sample tubes with prefilled sample extraction buffer
- 10 saliva collectors
- 1 tube stand
- 1 Instruction pamphlet

### Additionally required:

- 1 timer

### TEST PREPARATION

Let test cassette and test components stand at a room temperature (15°C to 27°C) before performing the test. Lay all the supplied materials on a clean, dry and flat surface.

### TEST PERFORMANCE


Read the instructions for use completely before performing the test. A step-by-step instruction is given on the next page and describes the test procedure.

### EVALUATION OF TEST RESULTS

To read the test results simply determine whether a line is present or absent at the Control (C) position. It does not matter how strong or weak a Control line (C) is.

### POSITIVE TEST RESULT


If a colored line is visible in the Test line area (T) and a colored line in the Control line (C) the test result is positive. The test result means that SARS-CoV-2 antigen are detectable in your saliva sample. The detection of these antigens indicates with a high probability of infection with the novel coronavirus.

\*Note: The thickness of the line is insignificant; any reddish color in the Test line (T) should be considered a positive result. The positive test result must be confirmed by PCR.

### NEGATIVE TEST RESULT

If only a colored line is visible in the Control line area (C) the test result is negative.



The test result indicates that there is no or too little SARS-CoV-2 Antigen in the saliva sample and at the current time there is probably no infection with the novel coronavirus.

**A negative result does not preclude SARS-CoV-2 infection.**

False negative results can be from incorrect sampling, incorrect execution of the test, or insufficient virus in the sample.

### INVALID TEST RESULT

If there is no Control line (C) or only a Test line (T) in the result window, the test did not run correctly and the results are not valid.

It is important that you carefully follow the instructions for the test. You should test again with a new saliva sample and a new test.

### PERFORMANCE EVALUATION

	Reference test: PCR		
	Positive	Negative	Total
NanoRepro SARS-CoV-2 Saliva Antigen Test	99	4	103
	6	396	402
	105	400	505

Sensitivity: 94,29% (87,98%-97,87%)\*  
 Specificity: 99,00% (97,46%-99,73%)\*  
 Trueness: 98,02% (96,39%-99,05%)\*  
 Accuracy: 96,12% (90,31%-98,50%)\*

\* 95 % Confidence interval

No cross-reactivity with the human coronaviruses (229E, OC43, NL63, SARS und MERS), Adenovirus, human Metapneumovirus, Parainfluenza virus 1,2,3 and 4, Influenza A, Influenza B, Enterovirus, Respiratory syncytial virus, Rhinovirus, Haemophilus influenzae, Streptococcus pneumoniae, Bordetella pertussis, Mycoplasma pneumoniae, Chlamydia pneumoniae, Legionella pneumophila, Candida albicans and pooled human nasal wash was observed when tested using the NanoRepro SARS-CoV-2 Saliva Antigen Test.

Following substances do not affect the performance of the NanoRepro SARS-CoV-2 Saliva Antigen Test: Whole Blood, Mucin, Nasal Drops (Sodium Chloride), Nasal Spray (Fluticasone Propionate, Gluconic Acid Zinc, Fluconazole, Oxymetazoline, Cromolyn), Sore Throat Phenol Spray, Throat Lozenge (Benzocaine, Menthol), Anti-viral Drug (Tamiflu, Ribavirin), Tobramycin 4.0

### WARNINGS AND IMPORTANT INFORMATION

- This kit is a qualitative detection, which cannot determine the exact content of antigen.
- The test is intended for use outside the body only.
- Not to be taken internally. Avoid sample buffer contact with skin and eyes.
- Keep out of the reach of children. Any child under age 18 shouldn't perform the test without parental guidance, or professional aid.
- Protect from sunlight, do not freeze. Store in a dry place between 2°C and 30°C. Do not use after the expiration date printed on the package.
- Not following the exact instructions can affect the outcome of the test. The final diagnosis must be confirmed by a physician.

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Follow instructions	In vitro diagnostic medical device (for external use)	Best before (See imprint on package)
Store at 2°C - 30°C Do not freeze	Contains sufficient for 10 tests	Do not reuse
Manufacturer	Keep dry	LOT Batch number (See imprint on package)
REF Reference number		

REF: B67000

Instructions English  
Revision from 2021-02 (Rev. 00)

NanoRepro AG  
Untergasse 8  
D-35037 Marburg  
www.nanorepro.com

## For healthcare professionals STEP-BY-STEP-INSTRUCTION

• Do not use the test if the packaging is damaged. Do not use broken test components.

• All test components are only intended to be used for this test. Do not reuse the test or test components.

• The test should be carried out immediately or within one hour after opening the foil pouch (20-30°C, humidity <60%).

• After human saliva sample is collected, the saliva should be processed as soon as possible and tested within 1 hour. If it cannot be tested immediately, it can be stored at 2-8°C for 72 hours.

• Poor vision, color blindness or poor lighting may affect your ability to interpret the test correctly.

• All specimens and reagents should be considered potentially hazardous and handled in the same manner as an infectious agent after use.

• A negative result does not rule out the infection of a SARS-CoV-2 infection. Therefore, the test should not be used as the only reference for the clinical diagnosis. The result must be confirmed by the PCR.

• The test is not validated on specimens from pregnant women.

### LIMIT OF DETECTION AND LITERATURE

The limit of detection of ZuhauseTest Corona has been evaluated at 1.62 x 10<sup>3</sup> TCID<sub>50</sub>/ml.

1.) Nanshan Chen\*, Min Zhou\*, Xuan Dong\*, Jieming Qu\*, Fengyun Gong, Yang Han, Yang Qiu, Jingli Wang, Ying Liu, Yuan Wei, Jia'an Xia, Ting Yu, Xinxin Zhang, Li Zhang Epidemiological and clinical characteristics of 99 cases of 2019 novel coronavirus pneumonia in Wuhan, China: a descriptive study. LANCET. January 29, 2020.

2.) World Health Organization (Coronavirus disease 2019) [https://www.who.int/emergencies/diseases/novel-coronavirus-2019/technical-guidance/naming-the-coronavirus-disease-\(covid-2019\)-and-the-virus-that-causes-it](https://www.who.int/emergencies/diseases/novel-coronavirus-2019/technical-guidance/naming-the-coronavirus-disease-(covid-2019)-and-the-virus-that-causes-it) (Zugriff am 27.03.2020)

3.) World Health Organization (Coronavirus disease 2019) [https://www.who.int/health-topics/coronavirus#tab=tab\\_1](https://www.who.int/health-topics/coronavirus#tab=tab_1) (Zugriff am 27.03.2020)


4.) Saliva as a Candidate for COVID-19 Diagnostic Testing: A Meta-Analysis Laszlo Márk Czumbel, Szabolcs Kiss, Nelli Farkas, Ivan Mardel, Anita Hegyi, Akos Nagy, Zsolt Lohinai, Zsolt Szakacs, Peter Hegyi, Martin C. Steward and Gabor Varga

### STEP 1


Do not eat, drink, smoke or chew gum for 30 minutes before collecting saliva sample. This may affect the test results.

### STEP 4


Place the tip of tongue against the mucous membrane of the oral cavity to enrich saliva, and then swallow saliva to the throat without ingesting saliva to esophagus, followed by returning saliva to the mouth. Repeat the process for 5 times, then spit saliva into a saliva collector. Avoid bubbles while collecting saliva.

### STEP 7


Close the cap of the sample tube.

### STEP 2


Open the sealed pouch and remove the test cassette. Lay it face up on a clean, dry and flat surface.

### STEP 5

Vertically hold the saliva collector tube and unscrew the saliva collector funnel to remove it from the tube. Discard the funnel. Take sample tube with prefilled sample extraction buffer and remove the white cap of the sample tube.

### STEP 8


Add 3 full drops of the mixed solution vertically into the sample well (S) of the test cassette.

Please note that there should be no liquid applied to the result windows marked with the letters (T) and (C). Do not touch or move the test cassette after adding the drops to the sample well (S).

Read the result 15-20 minutes after adding the sample. Result got after 20 minutes is invalid.


### STEP 3

Screw saliva collection funnel on to the saliva collector tube.

Attention: This tube does not contain liquid!

### STEP 6


Draw 2/3 dropper range of saliva sample from the saliva collector tube with the help of pipette and transfer the entire contents into the extraction buffer of the sample tube.

To thoroughly mix the saliva sample with extraction buffer, draw the solution up and down with the pipette about 10 times.


Attention: Please avoid bubbles!

### RESULTS



#### POSITIVE

#### NEGATIVE

#### INVALID

Rapid test for the detection of SARS-CoV-2 Antigen