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EINFÜHRUNG

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



INSTRUCTIONS FOR USE

INTRODUCTION

FertiQUICK is a rapid test for use at home to detect the concentration of sperm in semen. The self-test shows whether the sperm count is considered within normal limits. Other fertility factors are not detected by the test.

As an immunochromatographic rapid test FertiQUICK detects a component of the sperm head (ACRV1), which is then bound by colloidal gold-labelled antibodies. This complex binds to the test line (T) and forms a visible red line. The test includes a built-in control system realized by a control line (C). The C-line confirms sufficient amount of sample volume as well as correct performance of the test.

An explanation of how to read and interpret the test results is given in the "Evaluation of test result" section. It is important to fully understand the entire instruction pamphlet before performing the test.

A **step-by-step instruction** is given on the **next page** and describes the test procedure.

TEST CONTENTS

- 1 test cassette in a sealed pouch (SP-10/ 15 Mio./ml)
- 1 solution bottle
- 1 semen transfer syringe

MedNet GmbH, Borkstraße 10, 48163 Münster

- 1 semen collection cup
- 1 instruction pamphlet

Additionally required:

- 1 timer

TEST PREPARATION

Please note the following before you obtain the semen sample by manual stimulation (masturbation):

1. Wait at least 48 hours, but not more than 7 days after the last ejaculation.
2. Do not use any lubricants or lotions.
3. Do not use a condom.
4. Place all of the test components and a timer on a flat surface. Remove the semen transfer syringe **but not yet** the test cassette from the package.

TEST PERFORMANCE

Let the test cassette and the solution bottle stand at the room temperature (15°C to 27°C) before performing the test.

1. Collect the semen sample into the semen collection cup. If you do lose some of the semen, discard the sample. Rinse the cup **ONLY** with tap water. Wait at least 48 hours, but not more than 7 days (see 1 above.) and collect a fresh sample.
2. Gently mix the semen sample in the collection cup and let it stand for 20-60 min. If the collected semen sample has visibly liquefied, perform the test within 3 hours.
3. Gently stir the semen sample with the semen transfer syringe in the cup about 10 times until it is well mixed.
4. Insert the semen transfer syringe into the semen sample avoiding any solid or sticky material within the semen. Pull the plunger to draw the sample into the semen transfer syringe **until it**

reaches the 0,1ml line on the syringe. Avoid getting air bubbles in the semen transfer syringe. If this happens, push the semen back out completely and then draw semen into the semen transfer syringe again.

5. Open the seal of the solution bottle.
6. Insert the semen transfer syringe into the solution bottle and push the plunger to add the entire semen sample (0,1ml) to the solution bottle.
7. Screw the cap back onto the solution bottle and gently mix the contents thoroughly by turning the sample tube 5-10 times upside down.
8. Let the solution bottle stand for **2 minutes** to create the semen mixture. If this does not happen, turn the solution bottle upside down again. Open the foil pouch and remove the test cassette. Lay it face up on a clean, dry and flat surface.
9. Twist off the small cap on the tip of the solution bottle cap. Hold the bottle with the semen mixture straight up and down over the test cassette and squeeze gently to add **exactly 2 drops** to the sample well (S). **Please note, that there should be no liquid applied to the result windows marked with the letters (T) and (C).**

Read the result 5-10 minutes after adding the sample. Low values may be missed before the 5 minutes have elapsed, and false positive results may occur after more than 10 minutes.


EVALUATION OF TEST RESULTS

To read the test results simply determine whether a line is present or absent at the control (C) position. It does not matter how strong or weak a Control line (C) is.

POSITIVE TEST RESULT

If there are any line at all at the **Control (C) position** and the **Test (T) position**, the test result is positive. It does not matter how strong or weak a Test line (T) is. The positive result means that the sperm count is at least 15 million sperm per milliliter and the sperm count level is considered normal to father a child.

However, a positive test result alone does not prove the fertility. It is important to understand that several factors besides the sperm count could affect the ability to father a child. So if you and your partner have been trying to conceive a child without success for several months, you should both see a physician for the comprehensive fertility evaluation.


NEGATIVE TEST RESULT

If there is **only Control line (C)** and **NO Test line (T)**, the test result is negative. The sperm count is less than 15 million per milliliter. However, a negative test result alone does not prove the infertility.

Some men have sperm counts below this range, so they still are able to father a child naturally. It is also important to understand that the sperm count can vary from day to day. So it is possible that you might get a positive result if you were to wait a while and test again. However, we recommend to talk to physician about the test result and



to have a complete semen analysis. It is particularly important to consult a physician if you and your partner have been trying to conceive a child for a several months or more without success.

INVALID TEST RESULT

If there is **no Control line (C) or only a Test line (T)** in the result window, the test did not run correctly and the results are not valid.

It is important that you carefully followed the instructions for the test. You should test again with a new sample and a new test. Wait at least 48 hours, but not more than 7 days after your last ejaculation to retest.

PERFORMANCE EVALUATION

| | | Reference Test | | |
|-------------|--------------|----------------|------------|------------|
| | | Positive | Negative | Total |
| Ferti-QUICK | Positive | 130 | 4 | 134 |
| | Negative | 0 | 111 | 111 |
| | Total | 130 | 115 | 245 |

Sensitivity: 100.00% Specificity: 96.52%
Trueness: 98.37% Accuracy: 97.01%

WARNINGS AND IMPORTANT INFORMATION

- The test is intended for use outside the body only.
- Not to be taken internally. Avoid contact of sample dilution buffer with skin and eyes.
- Keep out of reach of children.
- Protect from direct sunlight, do not freeze. Store in a dry place between 2 °C and 30 °C.
- All test components are only intended to be used for this test. Do not reuse the test or test components!
- Do not use after the expiration date printed on the package.
- Not following the exact instructions can affect the outcome of the test. The definitive diagnosis must be confirmed by a physician.
- This test only assesses the sperm concentration and can not provide further information about fertility.
- Do not use the test if the packaging is damaged. Do not use broken test components.
- This test cannot be used as a method of birth control, paternity proof or as a protection against sexually transmitted diseases.
- The test should be performed immediately or within one hour after opening the sealed pouch.
- Poor vision, color blindness or poor lightening may affect your ability to interpret the test correctly.
- All test components can be disposed of in the household waste.

BIOLOGICAL REFERENCE RANGE AND LITERATURE

According to WHO-manual no. 5 the minimum amount of sperm required for successful reproduction is 15 million/mL.

1) WHO laboratory manual for the examination and processing of human semen - 5th ed. World Health Organization

2) M.A. Coppola1, K.L. Klotz2, K.-a. Kim3, etc. SpermCheck Fertility, an immunodiagnostic home test that detects normozoospermia and severe oligozoospermia.

www.fertiquick.de

Explanation of symbols:

| | | |
|------------------------------------|---|--|
| Follow instructions | In vitro diagnostic medical device (for external use) | Best before (See imprint on package) |
| Store at 2°C - 30°C Do not Freeze. | Content sufficient for 1 test | Do not reuse |
| Manufacturer | Sterilisation by ethylene oxide | Lot Batchnumber (See imprint on package) |
| Catalogue number | Keep dry | European Representative |

Instructions English
Revision from 2021-09 (Rev. 00)

NanoRepro AG
Untergasse 8
D-35037 Marburg
www.nano.org

REF 820001


For the love of life.

For Self-testing

STEP-BY-STEP INSTRUCTION


For Self-testing

STEP 1


Collect the semen sample in the collection cup provided. Gently mix the semen sample in the collection cup.

STEP 2

Let semen stand for **20-60 minutes** at the room temperature (15°C-30°C), until the sample has visibly liquefied.


STEP 3


Use semen transfer syringe to gently stir the semen sample 10 times.

TEST RESULT



POSITIVE TEST RESULT


NEGATIVE TEST RESULT

INVALID TEST RESULT


STEP 4


Pull plunger to draw liquefied semen **exactly up to the 0,1ml** line on the syringe.

STEP 5


Open the solution bottle by removing the seal of the buffer tube.

STEP 6


Insert the semen transfer syringe into the solution bottle and push the plunger gently to add all of the semen to the solution.

STEP 7


Screw the cap back onto the solution bottle and gently mix the contents thoroughly by turning the sample tube upside down (5-10 times). **Do not shake!**

STEP 8

Let the solution bottle containing the semen mixture stand for 2 minutes.

STEP 9

Open the foil pouch containing the test cassette and lay it face up. Twist off the small cap on the tip of the solution bottle cap and squeeze gently to add **exactly 2 drops** of the semen mixture to the sample well (S) of the test cassette.

Read the result **5-10 minutes** after adding the sample. **Low values may be missed before the 5 minutes have elapsed, and false positive results may occur after more than 10 minutes.**

www.fertiquick.de

For the love of life.