GEBRAUCH

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
INSTRUCTIONS FOR USE

INTENDED USE

GraviQUICK® is a self-performing immunoassay designed for the qualitative determination of human chorionic gonadotropin (hCG) in urine. **GraviQUICK®** is designed for both professional and lay users.

SUMMARY AND EXPLANATION

Human Chorionic gonadotropin (hCG) is a glycoprotein hormone secreted by the developing placenta shortly after fertilization. During normal pregnancy, hCG can be detected in urine as early as 7 days following conception, doubling every 1.3 to 2 days. At the time of the last missed menstrual period, urine hCG levels are about 100 mIU/ml with peak levels of 100,000 to 200,000 mIU/ml at the end of the first trimester. The presence of hCG soon after conception and its subsequent increase in concentration during early gestational growth make it an ideal marker for the early detection of pregnancy.

PRINCIPLE

GraviQUICK® is a rapid one-step test, based on an immunochromatographic technology. A membrane with an absorbent pad overlaps a strip of fiber glass paper that is impregnated with a lyophilized colloidal conjugate of gold particles and monoclonal solid phase antibodies to hCG. Another absorbent pad at the end of the assay absorbs excess sample fluid. The urine sample is introduced into the device, and proceeds through the absorbent pad, then laterally onto a chromatographic membrane. As it contacts the membrane, the sample dissolves the lyophilized conjugate. In a reactive sample, the hCG antigen will bind to the antibodies in the colloidal solution. As the conjugate moves forward on the membrane, anti-hCG monoclonal antibody affixed on the test zone ("T") will bind the hCG-gold conjugate complex, forming a pink line ("T"). All samples will cause a pink colored line in the control zone ("C"). This line is formed by the binding of the polyclonal antibodies (Anti-mouse IgG) affixed onto the control zone to the sample-colloidal gold conjugate. Presence of this line indicates sample amount is appropriate and the test has been carried out correctly. In less than 5 minutes, levels of hCG as low as 10 mIU/ml can be detected.

REAGENTS

Test device comprised colloidal gold coated with anti β-hCG monoclonal antibody; NC membrane coated with anti α-hCG monoclonal antibody and rabbit anti mouse IgG polyclonal antibodies.

WARNING

- Read directions for use carefully before performing this test. Pay attention to the position of the C and T line.
- Do not use beyond the labeled expiration date.
- Do not reuse the test devices. Discard it in the dustbin after single use.
- Do not use if pouch is damaged or opened.
- Do not touch the membrane on the strip.
- Once open the pouch, the test device should be used immediately. Prolonged exposure to ambient humidity will cause product deterioration.
- Treat urine samples and used devices as if they are potentially infectious. Avoid contact with skin.
- Ensure there is urine container before test.

TEST CONTENTS

- 1 test cassette and 1 pipette in foil pouch with 1 pack of desiccant (hCG 10mIU/ml)
- 1 instruction pamphlet

ADDITIONALLY REQUIRED

- 1 timer
- 1 urine cup
- 1 quality control material

Commercial control material or fresh samples from known pregnant and non-pregnant patients may be used with each lot as controls to verify the performance of the test. We recommend using control materials with hCG concentrations known to be near (above and below) the clinical cutoff.

STORAGE AND STABILITY

The test kit can be stored at normal temperature (4 – 30 °C, or 39 – 86° F) in the sealed pouch to the date of expiration. The test kits should be kept away from direct sunlight, moisture and heat. Do not freeze.

ASSAY PROCEDURE

Determination of test date


The test could be used from the 10th day after conception. In most cases the test is accurate at this date. However, because of variation in hCG concentration among women, or potential miscalculation of the last period, if there is a negative result, and pregnancy is still suspected, you should retest after a few days.

Specimen collection and handling


First morning urine specimen is preferred due to its high concentration of hCG. However, urine specimens collected at any time of the day may be used. If specimens cannot be assayed immediately, they can be stored at 2 – 8°C for up to 72 hours prior to testing. Urine specimens exhibiting visible precipitates should be settled to obtain a clear specimen for testing.

Test procedure

1. Bring the test pouch and urine to room temperature. To begin testing, open the sealed pouch by tearing along the notch. Remove the test from the pouch and use it as soon as once opened.
2. Using the sample dropper, draw about 3 drops (approximately 120 – 150 µl) urine sample from the specimen cup into the pipette, and dispense it into the sample well on the cassette (see the picture).
3. Wait for colored bands to appear. Depending on the concentration of hCG in the test specimen, positive results may be observed in as short as 40 seconds. However, to confirm Negative results, the complete reaction time (5 minutes) is required. So read the result within 5 minutes.
4. Do not read results after 10 minutes.

INTERPRETATION OF RESULTS

POSITIVE:
Two distinct color bands appear – one on the control and one on the test region. This indicates that pregnancy has been detected. The color intensity of the test bands may vary during pregnancy.

NEGATIVE:

Only one color band appears on the control line, and no band appears on the test region. This means there is no pregnancy detected.

INVALID:

There are no visible bands at all, or if one red band appears only on the test region. If this happens, repeat a new test. If test still fails, please contact the distributor with the lot number.

QUALITY CONTROL

Testing of quality control material to verify the performance of the test is recommended for each new lot of test. In addition, users should follow their state and local regulations and guidelines regarding GLP requirements. In addition, a built-in procedural control is being used designated as control region (C) on the test device to indicate that sufficient sample volume was added.

PERFORMANCE CHARACTERISTICS

Sensitivity

GraviQUICK® will display positive results with specimens containing hCG at the level close to or greater than 10 mIU/ml.

Specificity

The following compounds exhibited no interference when dissolved in urine. Non-cross reacting homologous hormones:

hFSH	1.000 mIU/ml (WHO 2nd IS)
hLH	500 mIU/ml (WHO 2nd IRP)
hTSH	1000 µIU/ml (WHO 2nd IRP)

Accuracy

Comparison studies on the **GraviQUICK®** with a legally marketed device were performed in various clinical reference laboratory. Positive and negative results were compared and the correlation was >99.9%.

Clinical Accuracy of GraviQUICK®:

	Reference method (Positive)	Reference Method (Negative)	Total
GraviQUICK®			
Positive	349	0	349
Negative	0	671	671
Total	349	671	1020

Diagnostic Sensitivity (Positive agreement) = 349/349 = 100%

Diagnostic Specificity (Negative agreement) = 671/671 = 100%

Overall Agreement: 1020/1020 x 100% = 100%

Non-cross reacting compounds:

Acetaminophen	20 mg/dl	Tetracycline	20 mg/dl	Ketone	20 mg/dl
Acetosal	20 mg/dl	Hemoglobin	20 mg/dl	Bilirubin	20 mg/dl
Salicylic Acid	20 mg/dl	Albumin	1000 mg/dl	Theelol	20 mg/dl
Ascorbic Acid	20 mg/dl	Glucose	20 mg/dl	Pregnanedione	20 mg/dl
Caffeine	20 mg/dl	Gentisic	20 mg/dl	Thiophene	20 mg/dl
Ampicillin	20 mg/dl				

Repeatability and Reproducibility:

Three lots of **GraviQUICK®** were used and 10 strips of each lot were tested in HCG standard solution with different concentrations prepared at 0 mIU/ml, 10 mIU/ml, 25 mIU/ml, 50mIU/ml and 100 mIU/ml.

hCG concentration	Lot 1		Lot 2		Lot 3		Variability
	P	N	P	N	P	N	
0 mIU	0	10	0	10	0	10	0%
10 mIU	10	0	10	0	10	0	0%
25 mIU	10	0	10	0	10	0	0%
50mIU	10	0	10	0	10	0	0%
100 mIU	10	0	10	0	10	0	0%

P: positive; N: negative

All samples were positive at and beyond the cut-off level. No different results were observed.

LIMITATION OF THE PROCEDURE

- Alcohol may interfere with test result. It is not recommended using the test after drinking.
- Occasionally specimens containing less than 10 mIU/ml for urine also yield positive results.
- A very early pregnancy containing an extremely low concentration of hCG can give a negative result. In this case, another specimen should be obtained at least 48 hours later and tested.
- hCG levels may remain detectable for several weeks after normal delivery, delivery by caesarean section, spontaneous abortion or therapeutic abortion.
- In cases where very high levels of hCG are present (> 500,000 mIU/ml) a false negative result can occur due to a "Prozone" effect. If pregnancy is still suspected, simply dilute specimen 1:1 with deionized water and retest.
- If a urine sample is too dilute (i.e.: low specific gravity) it may not contain a representative level of hCG. If pregnancy is still suspected, a first morning urine sample should be obtained and retested 48 hours later.
- As is true with any diagnostic procedure, the user should evaluate data obtained by the use of this kit in light of other clinical information and consult to the physicians for the final diagnosis of pregnancy before make any decision of medical relevance.

REFERENCE

1. Chard T. Pregnancy tests: a review. Hum Reprod. 1992 May; 7(5):701-10. Review.
2. American Pregnancy Association. Human Chorionic Gonadotropin (hCG): The Pregnancy Hormone.
3. Womens health Research. Dept. of OB/GNY, University of New Mexico. HCG Reference Service.

For further questions and additional literature please contact the manufacturer.

Explanation of symbols:

Follow instructions	In vitro diagnostic medical device (for external use)	Best before (See imprint on package)	Manufacturer	Keep away from sunlight	Caution (observe warnings)	Keep dry
Store at 4-30°C. Do not freeze!	Content sufficient for 1 test..	Do not reuse	Lot Batch number (see imprint package)	Do not use if package is damaged	Catalogue number	

CE 0483

Instructions English
Revision from April 2021 (Rev.02)

REF 200500

NanoRepro AG
Untergasse 8
D-35037 Marburg
www.nano.ag

For the love of life.