COVID-19/FLU A+B/ RSV/Adenovirus Antigen Rapid Test

"Test it quickly"

self-test for use at home

Fast
Reliable
Accurate

INTENDED USE

The test is a rapid chromatographic immunoassay for the qualitative detection of SARS-CoV-2, Influenza A virus, Influenza B virus, Respiratory Syncytial Virus (RSV) and Adenovirus antigens in human nasal swab specimens as an aid in the diagnosis of SARS-CoV-2, Influenza A, Influenza B, RSV and Adenovirus infections.
The test is intended for individuals who are suspected of being infected with SARS-CoV-2, Influenza A/B, RSV and/or Adenovirus.
The test is not automated. For self-testing use. For *in vitro* diagnostic use only.

SUMMARY

The novel coronaviruses belong to the β genus. COVID-19 is an acute respiratory infectious disease. People are generally susceptible. Currently, the patients infected with the novel coronavirus are the main source of infection; asymptomatic infected people can also be a source of infection. Based on the current epidemiological investigation, the incubation period is 1 to 14 days, mostly 3 to 7 days. The main manifestations include fever, fatigue and dry cough. Nasal congestion, runny nose, sore throat, myalgia and diarrhea are found in a few cases.[1]
Influenza (commonly known as 'flu') is a highly contagious, acute viral infection of the respiratory tract. Virus transmission occurs when a susceptible individual comes into contact with aerosols or respiratory fomites from an infected individual.[2] Influenza outbreaks occur each year during the fall and winter months. Type A viruses are typically more prevalent than type B viruses and are associated with most serious influenza epidemics, while type B infections are usually milder.
Respiratory Syncytial Virus (RSV), which causes infection of the lungs and breathing passages, is a major cause of respiratory illness in young children. In adults, it may only produce symptoms of a common cold, such as a stuffy or runny nose, sore throat, mild headache, cough, fever, and a general feeling of being ill. Most children with RSV infection, both those who were hospitalized and those who were treated as outpatients, had no coexisting medical conditions or characteristics that significantly identified them as being at greater risk for severe RSV disease, except for being under 2 years of age.[3]
Human Adenoviruses comprise an important group of etiologic agents that are responsible for various diseases in adults and children, such as respiratory, ocular, gastroenteric, and urinary infections. In immunocompromised and organ-transplanted individuals, these agents can cause generalized infections.[4]

HOW DOES IT WORK

The test is a qualitative, lateral flow immunoassay for the detection of SARS-CoV-2, Influenza A and Influenza B viruses, RSV, and Adenovirus antigens in human nasal swab specimen. When the extracted specimen is added to the specimen well of the test cassette, the SARS-CoV-2, Influenza A virus, Influenza B virus, RSV, or Adenovirus antigens, if present in the specimen, react specifically with the virus antibodies coated onto the particles, forming a mixture. The mixture migrates up the membrane and reacts with the virus antibodies immobilized on the membrane, resulting in the generation of one or two colored lines in the test line region (T/A/B) of each test window. The presence of the colored line(s) in the test line region(s) indicates a positive result, while its absence indicates a negative result. To serve as a built-in quality control, a colored line will always appear in the control line region(C), indicating that the test cassette functions well.


REAGENTS

The test contains Mouse IgG coated colloidal gold particle, Rabbit IgG coated colloidal gold particle, SARS-CoV-2 antibody coated latex particle, Influenza A antibody coated colloidal gold particle, Influenza B antibody coated colloidal gold particle, RSV antibody coated colloidal gold particle, Adenovirus antibody coated colloidal gold particle, and contains Goat anti-Mouse IgG antibody, Goat anti-Rabbit IgG antibody, SARS-CoV-2 antibody, Influenza A antibody, Influenza B antibody, RSV antibody, Adenovirus antibody coated on the membrane.

KIT COMPONENTS

Materials Provided

Components	**Kit size**	**1 T/kit**	**5 T/kit**
	Test cassette(s)	**1**	**5**
	Package insert	**1**	**1**
	Sterile swab(s)	**1**	**5**
	Buffer 0.55mL (0.02% Proclin 300)	**1**	**5**
	Tube holder	**On the box**	

Materials Required But Not Provided

- Timer

SPECIMEN COLLECTION

Please wash and dry your hands thoroughly. Remove the sterile swab from the pouch.
Do not touch the soft tip of the swab. Nasal swab specimen collection:

1. Insert the swab into one nostril (Approx. 2 cm deep).
2. Slowly rotate the swab for 5-10 times, rubbing it against the nasal wall.


Note:
This may feel uncomfortable. Do not insert the swab any deeper if you feel strong resistance or pain.
When the nasal mucosa is damaged or bleeding, nasal swab collection is not recommended.
If you are swabbing others, please wear a face mask. With children, you may not need to insert the swab so deep into the nostril. For very young children, you may need another person to steady the child's head while swabbing.
3. Remove the swab and repeat the collection in the other nostril using the same swab, avoiding excessive or highly viscous nasal discharge.

Note: Test should be performed immediately after the specimen has been collected.

HOW TO DO THE TEST


Allow the test cassette, the specimen and buffer to equilibrate to room temperature (15-30°C) prior to testing.

1

Remove the cover of the Buffer tube and insert the tube into the dotted circle on the packaging box or into the tube holder.

2


Place the swab specimen immediately into the tube with buffer.

3


Press the swab specimen head against the inside of the tube and rotate the swab for 10-15 seconds.

4


Remove the swab while squeezing the swab head against the inside of the tube to release as much liquid as possible from the swab.
Note 1: Discard the swab in compliance with local regulations.
Note 2: Test should be performed immediately after the specimen is treated with buffer.

5

Fit the tube tip onto the tube.

6

Take the test cassette out of the foil pouch and place it on a flat and level surface. Invert the tube and add **3 drops** of solution into each specimen well (S) respectively.
Note: If a drop of solution added contains air bubbles, add another drop to the well.

7

Start a timer and read the result **10 minutes** later after adding the specimen.
Do not interpret the result **after 20 minutes**.
Note: Do not touch or move the test cassette during the waiting time.

8


When finished, dispose of all used components according to local regulations.

READING THE RESULTS

POSITIVE
* **COVID-19/RSV/Adenovirus:**
Two colored lines appear in the COVID-19/RSV/Adeno window. One colored line should be in the control line region (C) and the other colored line should be in the test line region (T).

POSITIVE
* **FLU A:**
Two colored lines appear in the FLU A+B window. One colored line should be in the control line region (C) and the other colored line should be in the FLU A region (A).

POSITIVE
* **FLU B:**
Two colored lines appear in the FLU A+B window. One colored line should be in the control line region (C) and the other colored line should be in the FLU B region (B).

POSITIVE
* **FLU A and FLU B:**
Three colored lines appear in the FLU A+B window. One colored line should be in the control line region (C) and two colored lines should be in the FLU A region (A) and FLU B region (B) respectively.

NEGATIVE
Only one colored line appears in the control line region (C). No colored line appears in the test line regions (T/B/A). It indicates that there are no SARS-CoV-2/Influenza A/Influenza B/RSV/Adenovirus antigens in the specimen or that their amount is too small to be detected.

INVALID
The line in the control line region (C) fails to appear. Incorrect testing procedure or test reagent deterioration are the most likely reasons. Please review the testing procedure and repeat the test with a new test kit. If the problem persists, discontinue using the kit immediately and contact your local distributor.
***NOTE:** The intensity of the colored line in the test line region (T/B/A) varies based on the amount of SARS-CoV-2, Influenza A, Influenza B, RSV, Adenovirus antigen present in the specimen. Thus, any shade of a colored line in the test line region (T/B/A) should be considered a positive result.

STORAGE AND STABILITY

1. Store the kit at 2-30 °C. DO NOT FREEZE.
2. Keep the kit away from sunlight, moisture and heat.
3. Do not use the kit beyond the expiration date.
4. Please use the test cassette within **one hour** after removing it from the foil pouch.

WARNINGS AND PRECAUTIONS

Please read all the information in this package insert before performing the test.

1. This kit is for *in vitro* diagnostic use only.
2. Do not eat, drink or smoke in the area where the specimens or kits are handled.
3. Do not drink the buffer in the kit. Carefully handle the buffer and avoid it contacting skin or eyes. Rinse with plenty of running water immediately if contacts.
4. Do not use the kit if any of the components or packaging is damaged.
5. **Use the kit only once** and follow the test procedures strictly. Do not dismantle the test cassette or touch the test window of the test cassette.
6. Keep the kit out of the reach of children. Testing for children should be conducted by an adult.
7. Handle all specimens as if they contain infectious agents. The used tests and specimens may be infectious and should be discarded according to local regulations.
8. Wash hands thoroughly before and after handling.
9. Please ensure that an appropriate amount of specimen is used for testing. Too much or too little specimen may lead to deviation of results.
10. Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority.
11. The components provided in the kit are approved for use with this test. Do not use components from any other commercial kit.
12. This test provides only a preliminary result and should not be considered a substitute for medical advice. Always seek guidance from a qualified healthcare provider before making any health-related decisions.

LIMITATIONS

1. The **SPECIMEN COLLECTION**, **HOW TO DO THE TEST** and **READING THE RESULTS** must be followed closely. Failure to follow the procedures may lead to inaccurate test results.
2. Children aged 3–14 years should be tested by an adult. Do not conduct the test if you do not feel confident testing a child. If the child experiences pain, stop the test immediately. This kit is not suitable for children under 3 years old.
3. If the test result is negative and clinical symptoms persist, it is probably because the amount of the virus in very early infection stage is too small to be detected, and it is recommended to do the test again with a new kit or with a molecular diagnostic device to rule out infection.
4. Negative results do not rule out SARS-CoV-2 infection, particularly in those who have been in contact with the virus. Follow-up testing with a molecular diagnostic technique should be considered to rule out infection in these individuals.
5. Positive results of COVID-19 may be due to infection with non-SARS-CoV-2 coronavirus strains or other interference factors.
6. A positive result does not preclude an underlying co-infection with another pathogen, therefore the possibility of an underlying bacterial infection should be considered.
7. Excess blood or mucus on the swab specimen may interfere with test performance, and even may yield a false positive result.
8. A false negative result may be obtained if the concentration of the viruses present in the specimen is not adequate or below the detectable level of the kit, or if you fail to follow these procedures such as improper specimen collection or testing.
9. The kit will only indicate the presence or the absence of SARS-CoV-2, Influenza A, Influenza B, RSV and Adenovirus antigens in the nasal swab specimen. Neither the quantitative value nor the rate of increase in the concentration of these viruses can be determined by this qualitative test.

QUESTIONS & ANSWERS

1. How does the COVID-19/FLU A+B/RSV/Adenovirus Antigen Rapid Test work?
The kit can qualitatively detect SARS-CoV-2, Influenza A, Influenza B, RSV and Adenovirus antigens in nasal swab specimens through the specific antibodies it contains. A positive result indicates that SARS-CoV-2, Influenza A, Influenza B, RSV and/or Adenovirus is/are probably present in the specimen.
2. When should the kit be used?
You can do this test as an early diagnosis of SARS-CoV-2/Influenza A+B/RSV/Adenovirus infection when you (or the patient) have (has) the symptoms such as fever or chills, cough, fatigue, sore throat, headache, difficulty breathing, runny or stuffy nose, muscle and joint pain, wheezing, sneezing, decrease in appetite, new loss of taste or smell, nausea or vomiting, diarrhea, severe malaise (feeling unwell), bronchospasm. You can also use this test kit when you suspect that you (or the patient) may have been infected even you (or the patient) don't (doesn't) have any symptoms since there exists the chance of asymptomatic case.
3. Can the result be incorrect?
The result may be incorrect if inadequate or improper specimen is collected or the **COVID-19/FLU A+B/RSV/Adenovirus Antigen Rapid Test** gets wet before being used, or if the number of drops of solution added is less than 3 or more than 4. Besides, due to the immunological principle involved, there exist the chances of false results in rare cases. A consultation with the doctor is always recommended for such tests based on immunological principles.
4. How do I interpret the test results if the color and the intensity of the lines are different?
The color and the intensity of the lines have no importance for result interpretation. The lines should only be homogeneous and clearly visible. It should be considered as positive whatever the color or intensity of the test line is.
5. My test result is negative. Does that mean I'm not infected?
A negative result means that you (or the patient) are (is) probably not infected or that the viral load is too low to be recognized by the test.
In case of suspicion, repeat the test after 1-2 days, as the SARS-CoV-2/Influenza virus/RSV/Adenovirus cannot be precisely detected in all phases of an infection.
6. If the test result is positive, what should I do?
A positive result indicates that you (or the patient) probably have (has) a SARS-CoV-2, Influenza A, Influenza B, RSV or Adenovirus infection. Please see a doctor for medical aid.

PERFORMANCE CHARACTERISTICS

Accuracy
The test was evaluated with clinical specimens. **RT-PCR** was used as the reference method. Study results are presented in the tables below.

COVID-19 Test:

Method		RT-PCR		Total
		Positive	Negative	
Rapid Test	Positive	213	1	214
	Negative	5	1377	1382
Total		218	1378	1596
Sensitivity		97.7% (95% CI*: 94.7%~99.3%)		
Specificity		99.9% (95% CI*: 99.6%~100.0%)		
Accuracy		99.6% (95% CI*: 99.2%~99.9%)		

FLU A+B Test:

Method		FLU A			FLU B		
		RT-PCR		Total	RT-PCR		Total
		Positive	Negative		Positive	Negative	
Rapid Test	Positive	179	1	180	137	0	137
	Negative	9	939	948	8	983	991
Total		188	940	1128	145	983	1128
Sensitivity		95.2%(95% CI*:91.1%~97.8%)			94.5% (95% CI*: 89.4%~97.6%)		
Specificity		99.9%(95% CI*:99.4%~100.0%)			>99.9%(95% CI*:98.6%~100.0%)		
Accuracy		99.1% (95% CI*: 98.4%~99.6%)			99.3% (95% CI*: 98.6%~99.7%)		

RSV Test:

Method		RT-PCR		Total
		Positive	Negative	
Rapid Test	Positive	151	2	153
	Negative	6	833	839
Total		157	835	992
Sensitivity		96.2% (95% CI*: 91.9%~98.6%)		
Specificity		99.8% (95% CI*: 99.1%~100.0%)		
Accuracy		99.2% (95% CI*: 98.4%~99.7%)		

Adenovirus Test:

Method		RT-PCR		Total
		Positive	Negative	
Rapid Test	Positive	118	2	120
	Negative	3	789	792
Total		121	791	912
Sensitivity		97.5% (95% CI*: 92.9%~99.5%)		
Specificity		99.8% (95% CI*: 99.1%~100.0%)		
Accuracy		99.5% (95% CI*: 98.7%~99.8%)		

* Confidence Intervals

Sensitivity
The **COVID-19/FLU A+B/RSV/Adenovirus Antigen Rapid Test** can detect various virus strains of SARS-CoV-2, Influenza A/B viruses, RSV and Adenovirus. The confirmed LODs for each virus strain are presented in the table below.

Virus Strain	LOD
WHO SARS-CoV-2 antigen International Standard	30 IU/mL
BetaCoV/Wuhan/IPBCAMS-WH-01/2019	78 $TCID_{50}$/mL
A/Sydney/5/2021(H1N1)	50 $TCID_{50}$/mL
A/California/07/2009(H1N1)	5.25×10^3 $CEID_{50}$/mL
A/South Australia/69/2019(H3N2)	50 $TCID_{50}$/mL
A/Wisconsin/15/2009(H3N2)	6.2×10^4 $CEID_{50}$/mL
A/Anhui/1/2013(H7N9)	650 $TCID_{50}$/mL
A/Guangdong/17SF003/2016(H7N9)	4.5×10^3 $TCID_{50}$/mL
B/Austria/1359417/2021(Victoria)	50 $TCID_{50}$/mL
B/Florida/78/2015(Victoria)	1.4×10^3 $TCID_{50}$/mL
B/Darwin/58/2019(Yamagata)	100 $TCID_{50}$/mL
B/Brisbane/1/2007(Yamagata)	3.16×10^4 $TCID_{50}$/mL
A2(RSV type A)	2.5×10^3 $TCID_{50}$/mL
B WV/14617/85(RSV type B)	1.0×10^3 $TCID_{50}$/mL
GB(Adenovirus type 3)	1.0×10^3 $TCID_{50}$/mL
Gomen(Adenovirus type 7)	5.0×10^3 $TCID_{50}$/mL

The test can also detect strains of SARS-CoV-2, Influenza A/B viruses, RSV and Adenovirus as shown in the following table.

Virus Strains	Subtype	Concentration ($TCID_{50}$/mL)
SARS-CoV-2	Delta variants (B.1.617.2)	1.0×10^4
SARS-CoV-2	Delta variants (B.1.617.3)	1.0×10^4
SARS-CoV-2	Delta variants (AY Sub-lineages)	1.0×10^4
SARS-CoV-2	Omicron variants (B.1.1.529)	1.0×10^6
SARS-CoV-2	Omicron variants (BA.2)	1.0×10^6
A/Brisbane/02/2018	Influenza A(H1N1)	1.5×10^7
A/Victoria/2570/2019	Influenza A(H1N1)	2.2×10^5
A/Sydney/175/2022	Influenza A(H1N1)	2.9×10^6
A/Darwin/122/2018	Influenza A(H1N1)	3.3×10^5
A/Darwin/6/2018	Influenza A(H1N1)	1.1×10^4
A/Hangzhou/A158/2023	Influenza A(09H1N1)	5.25×10^6
A/Hangzhou/A419/2024	Influenza A(09H1N1)	2.63×10^6
A/Hubei/SWL1481/2023	Influenza A(09H1N1)	1.51×10^6
A/Brisbane/192/2017	Influenza A(H3N2)	3.1×10^6
A/Perth/9/2019	Influenza A(H3N2)	1.9×10^7
A/Darwin/9/2021	Influenza A(H3N2)	4.4×10^6
A/Victoria/6/2022	Influenza A(H3N2)	3.7×10^5
A/Sydney/1020/2018	Influenza A(H3N2)	6.4×10^7
A/Hubei/563/2019	Influenza A(H3N2)	5.9×10^6
A/Hangzhou/B870/2023	Influenza A(H3N2)	1.47×10^7
A/Hubei/VQ180801/2019	Influenza A(H5N1)	2.5×10^7
A/Anhui/VQ190501/2019	Influenza A(H5N1)	2.0×10^7
A/Hunan/VQ180901/2019	Influenza A(H5N1)	3.1×10^7
A/Hubei/c18259255/2018	Influenza A(H5N6)	1.0×10^5
A/LVS/2/2018	Influenza A(H5N6)	1.2×10^6
A/Hubei/c18255022/2018	Influenza A(H7N9)	1.0×10^5
A/Shanghai/02/2013	Influenza A(H7N9)	1.5×10^6
A/Hubei/c180510359/2018 (H9N2)	Influenza A (H9N2)	1.0×10^5
A/LVS/2/2018	Influenza A (H9N2)	1.4×10^6
B/Victoria/2110/19	Influenza B(Victoria)	1.7×10^6
B/Brisbane/35/2018	Influenza B(Victoria)	1.5×10^4
B/Victoria/705/2018	Influenza B(Victoria)	5.7×10^5
B/Hangzhou/B269/2024	Influenza B(Victoria)	6.67×10^6
B/Hubei/1832/2016	Influenza B(Victoria)	4.74×10^6
B/South Australia/67/2018	Influenza B(Yamagata)	2.4×10^4
B/Victoria/706/2018	Influenza B(Yamagata)	8.3×10^6
B/Hangzhou/B52/2018	Influenza B(Yamagata)	5.62×10^6
B/Hubei/24/2018	Influenza B(Yamagata)	2.16×10^6
Long	RSV type A	1.6×10^8
HRSV-A/Hep2/Hubei/2024-XN-335/2024	RSV type A	1.1×10^6
18537	RSV type B	1.6×10^7
RSVB/Human/Hubei/201710052/2018	RSV type B	1.0×10^6
Adenoid 71	Human adenovirus 1	1.0×10^7
Adenoid 6	Human adenovirus 2	1.0×10^6
RI-67	Human adenovirus 4	5×10^6
Clinical Virus Isolate	Human adenovirus 5	4×10^6
Tonsil 99	Human adenovirus 6	5×10^7
7A	Human adenovirus 7	3.16×10^5
Clinical Virus Isolate	Human adenovirus 11	3.1×10^6
Clinical Virus Isolate	Human adenovirus 14	3.1×10^5
Clinical Virus Isolate	Human adenovirus 55	1.6×10^6

Hook Effect
SARS-CoV-2 positive specimen of 1.0×10^7 $TCID_{50}$/mL, Influenza A virus positive specimen of 4.22×10^7 $TCID_{50}$/mL, Influenza B virus positive specimen of 1.78×10^7 $TCID_{50}$/mL, RSV positive specimen of 1.6×10^7 $TCID_{50}$/mL and Adenovirus positive specimen of 1.6×10^8 $TCID_{50}$/mL were tested and no high-dose hook effect was observed.

Precision-Repeatability
Precision-repeatability was determined by using specimens: negative, SARS-CoV-2 weak positive, SARS-CoV-2 strong positive, Influenza A weak positive, Influenza A strong positive, Influenza B weak positive, Influenza B strong positive, RSV weak positive, RSV strong positive, Adenovirus weak positive and Adenovirus strong positive specimens. The study was performed 5 replicates per day for 5 consecutive days by one operator using 1 lot of the test. No difference was detected in intra lot.

Precision-Reproducibility
Precision-reproducibility was determined using specimens: negative, SARS-CoV-2 weak positive, SARS-CoV-2 strong positive, Influenza A weak positive, Influenza A strong positive, Influenza B weak positive, Influenza B strong positive, RSV weak positive, RSV strong positive, Adenovirus weak positive and Adenovirus strong positive specimens. The study was performed 5 replicates per day for 5 consecutive days in 3 different sites using 3 lots of the test (one lot per site), and by three operators per site. No difference was detected between days, sites, lots and operators.

Cross-Reactivity
The test was evaluated with a total of 22 bacterial isolates at 1.0×10^7 CFU/mL or IFU/mL, and Pooled human nasal wash. None of the organisms listed below gave a positive result.

Chlamydophila pneumoniae	*Staphylococcus aureus subsp.aureus*
Mycoplasma pneumoniae	*Streptococcus pneumoniae*
Moraxella catarrhalis	*Streptococcus pyogenes*
Staphylococcus epidermidis	*Streptococcus salivarius*
Arcanobacterium	*Streptococcus sp Group F*
Corynebacterium	*Haemophilus influenzae*
Escherichia coli	*Bordetella pertussis*
Neisseria lactamica	*Legionella pneumophila*
Neisseria subflava	*Candida albicans*
Pseudomonas aeruginosa	*Klebsiella Pneumoniae*
Mycobacterium tuberculosis	*Pneumocystis jirovecii*
Pooled human nasal wash	

The test was evaluated with a total of 28 viral strains at concentrations listed in the table below. The study results showed that the test can specifically detect SARS-CoV-2, Influenza A/B viruses, RSV and Adenovirus with each strip respectively and there is no cross-reactivity with the virus strains below.

Virus Strain	Test Level
Adenovirus type 2	1.6×10^8 $TCID_{50}$/mL
Adenovirus type 3	1.6×10^8 $TCID_{50}$/mL
Adenovirus type 5	1.6×10^{10} $TCID_{50}$/mL
Adenovirus type 7	1.0×10^6 $TCID_{50}$/mL
Human coronavirus OC43	1.0×10^6 $TCID_{50}$/mL
Human coronavirus 229E	1.6×10^6 $TCID_{50}$/mL
Human coronavirus NL63	1.0×10^6 $TCID_{50}$/mL
Human coronavirus HKU1	1.0×10^6 $TCID_{50}$/mL
SARS-CoV-2	1.0×10^7 $TCID_{50}$/mL
MERS CoV Florida	1.05×10^6 $TCID_{50}$/mL
Human Rhinovirus 1A	2.2×10^7 PFU/mL
Human Rhinovirus 2	2.81×10^4 $TCID_{50}$/mL
Human Rhinovirus 14	1.58×10^6 $TCID_{50}$/mL
Human Rhinovirus 16	8.89×10^6 $TCID_{50}$/mL
Measles virus	8.9×10^5 $TCID_{50}$/mL
Mumps virus	2.8×10^6 $TCID_{50}$/mL
Parainfluenza virus 1	2.8×10^6 $TCID_{50}$/mL
Parainfluenza virus 2	1.6×10^6 $TCID_{50}$/mL
Parainfluenza virus 3	8.9×10^8 $TCID_{50}$/mL
Parainfluenza virus 4	2.8×10^6 $TCID_{50}$/mL
Respiratory syncytial virus	1.0×10^7 $TCID_{50}$/mL
Human Metapneumovirus	1.6×10^7 $TCID_{50}$/mL
Influenza A H1N1	1.3×10^7 $TCID_{50}$/mL
Influenza B	2.81×10^7 $TCID_{50}$/mL
Human Enterovirus 71	8.9×10^7 $TCID_{50}$/mL
Human herpesvirus 1	1.6×10^6 $TCID_{50}$/mL
Epstein-Barr virus	2.5×10^6 $TCID_{50}$/mL
Varicella-Zoster Virus	6.3×10^7 $TCID_{50}$/mL

Interfering Substance
The following substances, naturally present in respiratory specimens or that may be artificially introduced into the respiratory tract, were evaluated at specified concentrations. None of them were found to affect the test performance of the kit.

Substance	Concentration	Substance	Concentration
Whole Blood	40 µL/mL	Doxylamine Succinate	236 ng/mL
Mucin	5 mg/mL	Ephedrine	237 ng/mL
Budesonide Nasal Spray	200 µL/mL	Flunisolide	6.8 ng/mL
Zicam Allergy Relief Nasal Spray	200 µL/mL	Guaiacol glyceryl ether	3.5 ng/mL
Afrin Nasal Congestion Relief Pump Mist	200 µL/mL	Mupirocin	12 mg/mL
Amazon Basic Care Nasal Four Nasal Spray	200 µL/mL	Oxymetazoline	0.6 mg/mL
Albuterol	400 ng/mL	Phenylephrine	12 mg/mL
4-Acetamidophenol	200 µL/mL	Rebetol	4.5 µg/mL
Acetylsalicylic Acid	650 µg/mL	Relenza	282 ng/mL
Chlorpheniramine	145 ng/mL	Rimatadine	282 ng/mL
Dexamethasone	0.8 mg/mL	Seawater Nasal Spray	200 µL/mL
Dextromethorphan	1 µL/mL	Tobramycin	2.43 mg/mL
Diphenhydramine	5 µg/mL	Triamcinolone	1 mg/mL
Flonase Allergy Relief Nasal Spray	50 µL/mL	Phenol Glycerin Bacteriostasis Solution	150 µL/mL
Sodium Cromoglicate, Naphazoline hydrochloride and Chlorphenamine Maleate Nasal Spray	150 µL/mL	Chloraseptic Max Strength Sore Throat Lozenges	10 mg/mL
Oseltamivir Phosphate	5 mg/mL		

Summary of safety and performance is available on website:
https://ec.europa.eu/tools/eudamed/

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For the love of life.

Instructions English
Revision as of 2026-03-31(Rev.00)

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